

# CBJ Home Data Sheet



## Contents

<b>Introduction</b> .....	<b>3</b>
<b>Features and options</b> .....	<b>3</b>
<b>Usage</b> .....	<b>4</b>
Recommendations.....	4
Warnings.....	4
<b>Ordering Example</b> .....	<b>5</b>
<b>Functionality</b> .....	<b>6</b>
Buzzer indications.....	6
Learn mode function .....	6
<b>Front cover</b> .....	<b>7</b>
<b>Dimensions</b> .....	<b>8</b>
<b>System overview</b> .....	<b>9</b>
<b>Contacts</b> .....	<b>10</b>

## Introduction



The CBJ Home is a specially developed solution for patient lifts. The complete system consists of a control box and a battery enclosed in a single elegant module.

The system is approved according to medical safety standards and contains a series of features ensuring a safe comfortable lift, e.g. the CBJ Home is equipped with a soft-start function, electrical emergency lowering, emergency stop etc.

## Features and options

- Protection class IPX4
- Approx. 40 cycles per charge (as with CBJ)
- Mounting without bracket
- Soft start
- Colour grey RAL 7035
- Learn mode for adjustable current cut-off via a special hand control
- Emergency stop button
- Buzzer sounds if battery capacity is low (2-10 cycles left)
- Front cover with:
  - Green LED indication of handset activation
  - Yellow LED indication of charging
  - Button for electrical emergency lowering
- External charger possible CH01
- Internal charger possible (100 V AC – 240 V AC (50/60 Hz))

## Usage

CBJ Home with internal charger:	Nominal current draw max. 280 mA Power consumption (standby) max. 1.3 W Power consumption (charging) max. 12 W
Duty cycle:	Max. 10% or 2 min. continuous use then 18 min. without use
Operation temperature:	+5 °C to +40 °C
Storage temperature:	-10 °C to +50 °C
Relative humidity:	20% to 80% – non-condensing
Atmospheric pressure:	700 to 1060 hPa
Meters above sea level:	Max. 3000 meters
Approvals:	IEC60601-1 ANSI/AAMI ES60601-1

For details concerning CBJ Home compatibility, see Sales Backup chapter 3



## Recommendations


- Use only original LINAK mains cables to ensure proper connection to internal charger
- Always use fully charged batteries for learning mode procedures



## Warnings

- In order to avoid injury, the emergency-stop should be activated in transport situations.
- When “learn mode” is used, and channel 2 is pressed instead of channel 1, the CBJ Home will learn a new current limit of nearly 0 Amp. This will make it impossible to run the actuator with channel 2 until a new learn mode has been programmed.

## Ordering Example



## Functionality



### Front cover on the JUMBO Home™ with:

LED indication of charging (yellow) always lit when charging with an external charger, e.g. CH01. When charging with internal charger the LED will blink during charging and change to a constant yellow light when fully charged.

LED indication of handset activation (green) is also lit when connected to mains.

Button for electrical emergency lowering.

Battery must be charged, in order to operate via this button.

## Buzzer indications


- Battery: If the battery is low the buzzer will sound constantly when the hand control is activated.
- Overload: (Not for external charger versions produced before February 2010) CBJ Home will beep twice if the actuator is overloaded (over the current limit) - current will be cut-off and the actuator will stop. The actuator can be started again when the overload is removed.
- Overheated: (Not for external charger versions produced before February 2010) CBJ Home will beep twice if overheated - current will be cut-off and the actuator will stop. The actuator can be started again when the temperature has returned to normal.
- Learn mode: Single beep when current cut-off value is stored.

## Learn mode function

The patient hoist norm: EN10535 states that a patient hoist is not allowed to lift more than 1.5 times the max. weight it is intended for. For example, if the lift is intended for max. 100kg it should not be able to lift more than 150kg.

By using the learn mode function it is possible to programme the CBJ Home so that the control box automatically stops whenever the max. intended weight is exceeded.


## Front cover



Front cover nr. 00HBW073079  
to be fitted to hand control HB7x161-00


## Dimensions

### With external charger



Drawing number: 0899001


### With internal charger



Drawing number: 0899002



# System overview



# Contacts

## FACTORIES

### DENMARK - HEADQUARTERS

LINAK A/S  
PHONE: +45 73 15 15 15  
FAX: +45 74 45 80 48  
FAX (SALES): +45 73 15 16 13  
WWW.LINAK.COM

### CHINA

LINAK (SHENZHEN) ACTUATOR SYSTEMS, LTD.  
PHONE: +86 755 8610 6656  
PHONE: +86 755 8610 6990  
WWW.LINAK.CN

### SLOVAKIA

LINAK SLOVAKIA S.R.O.  
PHONE: +421 51 7563 444  
WWW.LINAK.SK

### THAILAND

LINAK APAC LTD.  
PHONE: +66 33 265 400  
WWW.LINAK.COM

### USA

LINAK U.S. INC.  
AMERICAS HEADQUARTERS  
PHONE: +1 502 253 5595  
FAX: +1 502 253 5596  
WWW.LINAK-US.COM  
WWW.LINAK-LATINAMERICA.COM

## SUBSIDIARIES

### Australia

LINAK Australia Pty. Ltd  
PHONE: +61 3 8796 9777  
FAX: +61 3 8796 9778  
E-mail: sales@linak.com.au  
www.linak.com.au

### Austria

LINAK Zweigniederlassung - Österreich (Wien)  
PHONE: +43 (1) 890 7446  
FAX: +43 (1) 890 744615  
E-mail: info@linak.de  
www.linak.at - www.linak.hu

### Belgium

LINAK Actuator-Systems NV/SA  
(Belgium & Luxembourg)  
PHONE: +32 (0)9 230 01 09  
E-mail: beinfo@linak.be  
www.linak.be - www.fr.linak.be

### Brazil

LINAK Do Brasil Comércio De Atuadores Ltda.  
PHONE: +55 (11) 2832 7070  
FAX: +55 (11) 2832 7060  
E-mail: info@linak.com.br  
www.linak.com.br

### Canada

LINAK Canada Inc.  
PHONE: +1 502 253 5595  
FAX: +1 416 255 7720  
E-mail: info@linak.ca  
www.linak-us.com

### Czech Republic

LINAK C&S s.r.o.  
PHONE: +42 058 174 1814  
FAX: +42 058 170 2452  
E-mail: info@linak.cz  
www.linak.cz - www.linak.sk

### Denmark - International

LINAK International  
PHONE: +45 73 15 15 15  
E-mail: info@linak.com  
www.linak.com

### Denmark - Sales

LINAK DANMARK A/S  
PHONE: +45 86 80 36 11  
FAX: +45 86 82 90 51  
E-mail: linak@linak-silkeborg.dk  
www.linak.dk

### Finland

LINAK OY  
PHONE: +358 10 841 8700  
E-mail: linak@linak.fi  
www.linak.fi

### France

LINAK FRANCE E.U.R.L.  
PHONE: +33 (0) 2 41 36 34 34  
FAX: +33 (0) 2 41 36 35 00  
E-mail: linak@linak.fr  
www.linak.fr

### Germany

LINAK GmbH  
PHONE: +49 6043 9655 0  
FAX: +49 6043 9655 60  
E-mail: info@linak.de  
www.linak.de

### India

LINAK A/S India Liaison Office  
PHONE: +91 120 4531797  
FAX: +91 120 4786428  
E-mail: info@linak.in  
www.linak.in

### Ireland

LINAK UK Limited (Ireland)  
PHONE: +44 (0)121 544 2211  
FAX: +44 (0)121 544 2552  
+44 (0)796 855 1606 (UK Mobile)  
+35 387 634 6554  
(Republic of Ireland Mobile)  
E-mail: sales@linak.co.uk  
www.linak.co.uk

### Italy

LINAK ITALIA S.r.l.  
PHONE: +39 02 48 46 33 66  
FAX: +39 02 48 46 82 52  
E-mail: info@linak.it  
www.linak.it

### Japan

LINAK K.K.  
PHONE: 81-45-533-0802  
FAX: 81-45-533-0803  
E-mail: linak@linak.jp  
www.linak.jp

### Malaysia

LINAK Actuators Sdn. Bhd.  
PHONE: +60 4 210 6500  
FAX: +60 4 226 8901  
E-mail: info@linak-asia.com  
www.linak.my

### Netherlands

LINAK Actuator-Systems B.V.  
PHONE: +31 76 5 42 44 40 /  
+31 76 200 11 10  
E-mail: info@linak.nl  
www.linak.nl

### New Zealand

LINAK New Zealand Ltd  
PHONE: +64 9580 2071  
FAX: +64 9580 2072  
E-mail: nzsales@linak.com.au  
www.linak.com.au

### Norway

LINAK Norge AS  
PHONE: +47 32 82 90 90  
E-mail: info@linak.no  
www.linak.no

### Poland

LINAK Polska  
LINAK Danmark A/S (Spółka Akcyjna)  
PHONE: +48 22 295 09 70 /  
+48 22 295 09 71  
E-mail: info@linak.pl  
www.linak.pl

### Republic of Korea

LINAK Korea Ltd.  
PHONE: +82 2 6231 1515  
FAX: +82 2 6231 1516  
E-mail: info@linak.kr  
www.linak.kr

### Slovakia

LINAK SLOVAKIA S.R.O.  
PHONE: +421 51 7563 444  
www.linak.sk

### Spain

LINAK Actuadores, S.L.U.  
PHONE: +34 93 588 27 77  
FAX: +34 93 588 27 85  
E-mail: esma@linak.es  
www.linak.es

### Sweden

LINAK Scandinavia AB  
PHONE: +46 8 732 20 00  
FAX: +46 8 732 20 50  
E-mail: info@linak.se  
www.linak.se

### Switzerland

LINAK AG  
PHONE: +41 43 388 31 88  
FAX: +41 43 388 31 87  
E-mail: info@linak.ch  
www.linak.ch - www.fr.linak.ch  
www.it.linak.ch

### Taiwan

LINAK (Shenzhen) Actuator systems Ltd.  
Taiwan Representative office  
PHONE: +886 2 272 90068  
FAX: +886 2 272 90096  
E-mail: sales@linak.com.tw  
www.linak.com.tw

### Turkey

LINAK İth. İhr. San. ve Tic. A.Ş.  
PHONE: +90 312 4726338  
FAX: +90 312 4726635  
E-mail: info@linak.com.tr  
www.linak.com.tr

### United Kingdom

LINAK UK Limited  
PHONE: +44 (0)121 544 2211  
FAX: +44 (0)121 544 2552  
E-mail: sales@linak.co.uk  
www.linak.co.uk

## DISTRIBUTORS

### Argentina

NOVOTEC ARGENTINA SRL  
PHONE: 011-4303-8989 / 8900  
FAX: 011-4032-0184  
E-mail: info@novotecargentina.com  
www.novotecargentina.com

### Colombia

MEM Ltda  
PHONE: +[57] (1) 334-7666  
FAX: +[57] (1) 282-1684  
E-mail: servicioalcliente@memltda.com.co  
www.mem.net.co

### India

Mechatronics Control Equipments India Pvt Ltd  
PHONE: +91-44-28558484, 85  
E-mail: bala@mechatronicscontrol.com  
www.mechatronicscontrol.com

### Indonesia

PT. HIMALAYA EVEREST JAYA  
PHONE: +6 221 544 8956,  
+6 221 544 8965  
FAX: +6 221 619 1925  
FAX (Sales): +6 221 619 4658  
E-mail: hejplastic-div@centrin.net.id  
www.hej.co.id

### Israel

NetivTech LTD  
PHONE: +972 55-2266-535  
FAX: +972 2-9900-560  
Email: info@NetivTech.com  
www.netivtech.com

### Singapore

Servo Dynamics Pte Ltd  
PHONE: +65 6844 0288  
FAX: +65 6844 0070  
E-mail: servodynamics@servo.com.sg

### South Africa

Industrial Specialised Applications CC  
PHONE: +27 011 466 0346  
E-mail: gartht@isagroup.co.za  
www.isaza.co.za

### United Arab Emirates

Mechatronics  
PHONE: +971 4 267 4311  
FAX: +971 4 267 4312  
E-mail: mechtron@emirates.net.ae

## Terms of use

LINAK® takes great care in providing accurate and up-to-date information on its products. However, the user is responsible for determining the suitability of LINAK products for a specific application. Due to continual development, LINAK products are subject to frequent modifications and changes. LINAK reserves the rights to conduct modifications, updates, and changes without any prior notice. For the same reason, LINAK cannot guarantee the correctness and actual status of imprinted information on its products.

LINAK uses its best efforts to fulfil orders. However, for the reasons mentioned above, LINAK cannot guarantee availability of any particular product at any given time. LINAK reserves the right to discontinue the sale of any product displayed on its website or listed in its catalogues or in other written material created and produced by LINAK, LINAK subsidiaries, or LINAK affiliates. All sales are subject to the 'Standard Terms of Sale and Delivery for LINAK A/S' available on LINAK websites. LINAK and the LINAK logotype are registered trademarks of LINAK A/S. All rights reserved.



WE IMPROVE YOUR LIFE